



Weight Calibration Service

Your UKAS Partner for Mass Units



Weight Calibration Service on a Range of Weights

If you are checking your weighing equipment on a regular basis with weights then you are already aware of the importance of accuracy of measurement. The same principles must therefore apply to the weights you are using to carry out these checks.

All weights are calibrated in a UKAS laboratory accredited to the international standard EN ISO IEC 17025 for testing laboratories.


Sartorius Mechatronics offers a recalibration service with a UKAS calibration certificate for weights ranging from:

1mg – 200g E2
1mg – 10Kg F1
1mg – 20Kg F2
1mg – 20Kg M1

Contact us now for a FREE quotation!

- Email: uk.customerservice@sartorius.com •
- Web: www.sartorius.co.uk •
- Web: www.sartorius.ie •
- Tel: 01372 737101 • Fax: 01372 722151 •
- Tel: 01 808 905 • Fax: 01 808 9388 •

We cover *All Brands* of Weights



How Do I Get My Weights Calibrated?

Simply contact Sartorius Mechatronics and tell us the amount and class of weights you require us to calibrate.

Sartorius will supply an estimate for the price of calibration, collection and return delivery.

Determining The Class of My Weights

E2 weights are constructed of stainless steel and have no markings on them.

F1 & F2 are constructed of stainless steel and will have the weight denomination stamped on the weight.

M1 are typically constructed of brass or cast Iron and have M stamped on them plus the denomination.

Sartorius Weight Calibration Contract

- Sartorius will provide a fixed cost per weight over a specified period
- Sartorius can provide a fixed cost of the pick-up and delivery of the weights
- A contract would be set-up and a planned calibration period agreed
- Sartorius will notify you when your weights are due for calibration and organise the complete process

Contact us now for a FREE quotation!

- Email: uk.customerservice@sartorius.com •
- Web: www.sartorius.co.uk •
- Web: www.sartorius.ie •
- Tel: 01372 737101 • Fax: 01372 722151 •
- Tel: 01 808 905 • Fax: 01 808 9388 •